

Pre-test Checklist- Verbal Instructions: AMCT

- Verify the examinee signed a consent form.
- If applicable, verify the examinee completed all pre-test documents.
- Inform the examinee the duration of the test is 30-35 minutes.
- Ensure the examinee does not require a restroom break.
- Ensure the examinee does not have any eye obstructions. If so, please remove. For example:
 - Baseball cap or hat with a rim that shadows the eyes
 - Hair that covers any part of the eyes
 - Thick, dark makeup
 - Long, fake eyelashes
 - Face mask positioned high on the cheeks
- Ask the examinee to turn off all electronic devices to avoid distractions such as notifications or vibrations.
- Briefly describe the test setup process: (1) finish verbal instructions, (2) load the test, (3) properly position the examinee, (4) calibrate eyes, and (5) run any optional diagnostic pre-tests.
- Explain that instructions and test questions will be read audibly; no text will appear onscreen. Mention to only click to proceed when voice tells indicates.
- Describe the testing process: (1) Pre-test instructions, (2) Practice test, and (3) Lie detection test. During the lie detection test, there are sequences of test questions followed by short breaks or rest periods.
- Explain all questions are TRUE/FALSE and question types include relevant topics and math.
- Describe how to answer with the mouse: TRUE (green button) and FALSE (red button). Physically show the examinee the red and green mouse buttons. Mention the mouse does not need to move to click the image of the mouse on screen. Mention the question must be completely read prior to clicking to answer.
- Explain that test questions must be answered quickly and accurately. It is ok to incorrectly answer or time out on a few questions. Forget mistakes and focus on the next question.
- Tell the examinee to blink normally.
- Ask the examinee to keep their chin in the chin rest for the entire test, except during the indicated breaks.
- Explain that breaks occur periodically and last 15 seconds; there is a countdown on the screen.
- Notify the examinee they will be monitored to ensure they are alert and cooperative. Both conditions are required to take the test. Voice will prompt them to refocus if their eyes close or get sleepy, which delays the test.
- Verify the examinee has no questions or concerns.
- Ask the examinee to sit.
- Do test set up process: load the test, properly position examinee, calibrate eyes (verbally explain), run any optional diagnostics (verbally explain), and start the test.

Sample pretest script:

This test will take about 30 to 35 minutes and should be uninterrupted. If interrupted, we may have to start over. If you need to visit the restroom, please do so now. Also, please power off any electronic devices to avoid distractions.

In this test, test instructions and questions will be read to you. No text will appear onscreen; instead, you should focus on a cluster of dots in the middle of the screen. When instructed by the voice, you will click to advance during the test. The test starts with pre-test instructions, then you will take a practice test to get used to answering the questions, and then you will take the lie detection test. In this test, you will be asked a series of true/false questions. Some questions will ask you about [insert 3 or 4 unique topics specific to your test. Example: illegal drug use, domestic violence, serious crimes, and terrorist associations]. And other questions will be basic math questions like $1 + 1 = 2$ which would be TRUE, or $2 + 3 = 7$, which would be FALSE. To answer true, please click the green mouse button. To answer false, click the red mouse button. You cannot answer until the voice has finished reading the entire question. You will know when you can answer when an image of a mouse appears onscreen. You do not need to move the cursor to click the onscreen mouse, just click the appropriate buttons on your physical mouse. It is important that you answer questions quickly and accurately. With that said, do not be worried if you incorrectly answer a few questions. That's normal. It does not mean you're going to fail the test. We do not expect perfection, but that you do your best. Also, please blink normally throughout the test.

During the test, please keep your chin in the chin rest, except during the specified break times. When there is a break, the test will prompt you to relax and then show a countdown onscreen. You have until the end of the countdown to relax, move around, or stretch your shoulders. After the countdown, place your chin in the chinrest to resume testing. I will monitor you during the test. It is important to remain alert and to cooperate with the test instructions. If your eyes wander or close, the screen will prompt you to refocus and present a 3-second countdown until you do, which could delay the test. Any questions?

Example of test set up script:

Before Calibration: *Please look at the dot and follow it as it moves along the screen.*

Before Pupil Light Reflex Diagnostic: *Please stare at the screen as it changes from light to dark.*

Before Digit Span Sequence Diagnostic: *You will be presented with a series of numbers. Memorize them. A keypad will appear on the screen, and you will use the mouse to enter in the numbers you memorized. You will gradually be presented with longer sets of numbers. Please do the best you can at accurately repeating them back.*